

# EIF4B antibody



Catalog Number: 110259

## Product name

EIF4B antibody

## Immunogen

[Human EIF4B Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EIF4B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 80 kDa

## Preparation

This antibody was obtained by immunization of EIF4B recombinant protein (Accession Number: NM\_001417). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110259(EIF4B antibody) at dilution of 1:300



IP Result of anti-EIF4B (IP:Catalog No:110259, 4ug; Detection:Catalog No:110259 1:300) with HeLa cells lysate 1080ug.