

Centrin 1 antibody



Catalog Number: 109259

Product name

Centrin 1 antibody

Immunogen

[Human CETN1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

Centrin 1 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in HL-60 cells, human heart tissue. Observed molecular weight by Western-blot: 20 kDa

Preparation

This antibody was obtained by immunization of Centrin 1 recombinant protein (Accession Number: NM_004066). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109259(CETN1 antibody) at dilution of 1:300



Immunofluorescent analysis of Hela cells, using CETN1 antibody Catalog No:109259 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).