

Product name

ERH antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ERH Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of ERH recombinant protein (Accession Number: NM_004450). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200

Validations

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110359(ERH antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:110359(ERH antibody) at dilution of 1:50 (under 10x lens)