

TREX1 antibody



Catalog Number: 116359

Product name

TREX1 antibody

Specificity

Human; other species not tested.

Antibody description

TREX1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of TREX1 recombinant protein (Accession Number: NM_016381). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116359(TREX1



Antibody) at dilution of 1:300

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116359(TREX1 Antibody) at dilution of 1:300



IP Result of anti-TREX1 (IP:Catalog No:116359, 4ug; Detection:Catalog No:116359 1:300) with HeLa cells lysate 1200ug.