

Anti-ASH2L Rabbit antibody



Catalog Number: 167359

Product name

Anti-ASH2L Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ASH2L Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ASH2L

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000-1/10000

IHC: 1/20


ICC/IF: 1/50

FC: 1/200

Validations



Immunocytochemistry of ASH2L(green) in HeLa cells using ASH2L Rabbit pAb at dilution 1/50, and



Western blot detection of ASH2L in C6,3T3,HeLa cell lysates using ASH2L Rabbit pAb(1:1000 diluted).Predicted band size:69kDa.Observed band size:80kDa.