

FYTTD1 antibody



Catalog Number: 110759

Product name

Polyclonal

FYTTD1 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

FYTTD1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 50 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of FYTTD1 recombinant protein (Accession Number: NM_032288). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110759(FYTTD1 Antibody) at dilution of 1:300