

Anti-Fatty Acid Synthase (2F9) Mouse antibody



Catalog Number: 168469

Product name

Anti-Fatty Acid Synthase (2F9) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Fatty Acid Synthase (2F9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Fatty Acid Synthase protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, IP, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/400

Validations




Western blot detection of Fatty Acid Synthase in HeLa, C6, 3T3, CHO-K1 and 293 cell lysates using Fatty Acid Synthase mouse mAb(dilution 1:1000). Predicted band size: 273kDa. Observed



band size: 273kDa.

Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-Fatty Acid Synthase mouse mAb (dilution 1:400).



Immunoprecipitation analysis of CHO-K1 cell lysates using Fatty Acid Synthase mouse mAb.