

WRN antibody



Catalog Number: 116969

Product name

Polyclonal

WRN antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

WRN Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells. Observed molecular weight by Western-blot: 162 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000553). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116969(WRN Antibody) at dilution of 1:300