

# DNMT3A antibody



Catalog Number: 109969

## Product name

DNMT3A antibody

## Immunogen

[Human DNMT3A Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DNMT3A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, Jurkat cells, MCF7 cells, mouse testis tissue, NIH/3T3 cells, PC-3 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 120-130kd

## Preparation

This antibody was obtained by immunization of DNMT3A recombinant protein (Accession Number: XM\_017003526). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109969(DNMT3A antibody) at dilution of



1:1500

IP Result of anti-DNMT3A (IP:Catalog No:109969, 4ug; Detection:Catalog No:109969 1:1000) with HeLa cells lysate 2000ug.