

Anti-Visfatin (2E11) Mouse antibody



Catalog Number: 168269

Product name

Mouse IgG1

Anti-Visfatin (2E11) Mouse antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description

WB: 1/500 - 1/2000

Visfatin (2E11) Mouse monoclonal antibody

ELISA: 1/10000

Preparation

Antigen: Purified recombinant fragment of PBEF1 (aa338-479) expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of PBEF/NAMPT in Jurkat, A549, MCF7, C6 and Hela cell lysates using PBEF/NAMPT mouse mAb (1:5000 diluted). Predicted band size: 56 kDa. Observed band size: 52 kDa.