

GAB3-Specific antibody



Catalog Number: 110769

Product name

GAB3-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GAB3-Specific Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue. Positive WB detected in mouse brain tissue, HepG2 cells, MCF7 cells. Observed molecular weight by Western-blot: 60-66 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001081573). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:50-1:500


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110769(GAB3-Specific antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:110769(GAB3-Specific Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:110769(GAB3-Specific Antibody) at dilution of 1:200 (under 40x lens).