

Anti-FAM98B Rabbit antibody



Catalog Number: 166769

Product name

Anti-FAM98B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

FAM98B Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human FAM98B

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000-1/5000


IHC: 1/20-1/100

ICC/IF: 1/20


Validations



Immunocytochemistry of FAM98B(green) in HeLa cells using FAM98B Rabbit pAb at dilution 1/50, and DAPI(blue)



Immunohistochemistry of FAM98B in paraffin-embedded Human lung cancer tissue using



FAM98B Rabbit pAb at dilution 1/20

Western blot detection of FAM98B in K562 cell lysates using FAM98B Rabbit pAb(1:1000 diluted). Predicted band size:37kDa. Observed band size:37kDa.