

CCRL2 antibody



Catalog Number: 109079

Product name

Polyclonal

CCRL2 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse; other species not tested.

Applications

Antibody description

ELISA, WB

CCRL2 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 47kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CCRL2 recombinant protein (Accession Number: NM_003965). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109079(CCRL2 antibody) at dilution of 1:400