

# ESRRB antibody



Catalog Number: 110479

## Product name

ESRRB antibody

## Immunogen

[Human ESRRB Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ESRRB Rabbit Polyclonal antibody. Positive WB detected in PC-13 cells, A549 cells. Positive IP detected in A549 cells. Observed molecular weight by Western-blot: 45-49 kDa

## Preparation

This antibody was obtained by immunization of ESRRB recombinant protein (Accession Number: XM\_017021086). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



PC-13 cells were subjected to SDS PAGE followed by western blot with Catalog No:110479(ESRRB antibody) at dilution of 1:1000



IP Result of anti-ESRRB (IP:Catalog No:110479, 4ug; Detection:Catalog No:110479 1:300) with A549 cells lysate 1040ug.