

MARK2 antibody



Catalog Number: 112479

Product name

MARK2 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MARK2 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive IP detected in rat brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 78 kDa, 80 kDa

Preparation

This antibody was obtained by immunization of MARK2 recombinant protein (Accession Number: BC008771). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112479(MARK2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human brain using Catalog No:112479(MARK2 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:112479(MARK2

MARK2 antibody



Catalog Number: 112479

antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-MARK2 (IP:Catalog No:112479, 4ug; Detection:Catalog No:112479 1:1000) with rat brain tissue lysate 4000ug.