

PPAR gamma antibody



Catalog Number: 107479

Product name

PPAR gamma antibody

Specificity

Human; other species not tested.

Antibody description

PPAR gamma Mouse Monoclonal antibody.
Positive WB detected in A431 cells. Positive IHC detected in human placenta tissue, human colon tissue. Observed molecular weight by Western-blot: 50-60kd

Preparation

This antibody was obtained by immunization of PPAR gamma recombinant protein (Accession Number: XM_011533841). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:107479(PPARG antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107479(PPARG Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107479(PPARG Antibody) at dilution of 1:50

PPAR gamma antibody



Catalog Number: 107479

(under 40x lens)