

USE1 antibody



Catalog Number: 116579

Product name

USE1 antibody

Immunogen

[Human USE1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

USE1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of USE1 recombinant protein (Accession Number: NM_018467). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116579(USE1 Antibody) at dilution of 1:300