

MYO16 antibody



Catalog Number: 112979

Product name

MYO16 antibody

Immunogen

[Human MYO16 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

MYO16 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 200-240 kDa

Preparation

This antibody was obtained by immunization of MYO16 recombinant protein (Accession Number: NM_015011). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112979(MYO16 Antibody) at dilution of 1:1000