

Anti-Calretinin Rabbit antibody



Catalog Number: 168779

Product name

Anti-Calretinin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Calretinin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Calretinin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000


ICC/IF: 1/20-1/50

FC: 1/20

Validations



Western blot detection of Calretinin in Rat Brain lysates using Calretinin Rabbit mAb(1:1000 diluted). Predicted band size:32kDa.Observed band size:32kDa.



Immunocytochemistry of Calretinin(green) in HeLa cells using Calretinin Rabbit mAb at dilution 1/50, and DAPI(blue)