

TMEFF2 antibody



Catalog Number: 116089

Product name

TMEFF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMEFF2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human prostate hyperplasia tissue. Positive WB detected in SH-SY5Y cells, HeLa cells, human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 45-48 kDa

Preparation

This antibody was obtained by immunization of TMEFF2 recombinant protein (Accession Number: NM_016192). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:116089(TMEFF2 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:116089(TMEFF2 Antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog

TMEFF2 antibody



Catalog Number: 116089

No:116089(TMEFF2 Antibody) at dilution of 1:100 (under 40x lens)