

# MUC16,CA125 antibody



Catalog Number: 112889

## Product name

MUC16,CA125 antibody

## Specificity

Human; other species not tested.

## Antibody description

MUC16,CA125 Rabbit Polyclonal antibody.  
Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 90-100 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_024690).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP, IHC

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:50-1:500

## Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112889(MUC16,CA125 antibody) at dilution of

1:300



IP Result of anti-MUC16,CA125 (IP:Catalog No:112889, 4ug; Detection:Catalog No:112889 1:600) with HeLa cells lysate 1200ug.



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using Catalog No:112889(MUC16,CA125 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen


# MUC16,CA125 antibody



Catalog Number: 112889

---

retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using Catalog No:112889(MUC16,CA125 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).