

FAM46B antibody



Catalog Number: 110589

Product name

FAM46B antibody

Immunogen

[Human FAM46B Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

FAM46B Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of FAM46B recombinant protein (Accession Number: NM_052943). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110589 (FAM46B Antibody) at dilution of 1:1000