

Anti-TUBGCP5 antibody



Catalog Number: 175989

Product name

Polyclonal

Anti-TUBGCP5 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat

Applications

Antibody description

WB, IHC-P

Rabbit polyclonal antibody to TUBGCP5

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tubgcp5. aa range:741-790

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1:; Western Blot analysis of various cells using GCP5 Polyclonal Antibody diluted at 1: 2000