

# XKR6 antibody



Catalog Number: 116989

---

## Product name

XKR6 antibody

## Specificity

Human; other species not tested.

## Antibody description

XKR6 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, Y79 cells. Observed molecular weight by Western-blot: 41-45 kDa

## Preparation

This antibody was obtained by immunization of XKR6 recombinant protein (Accession Number: BC024146). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116989(XKR6 Antibody) at dilution of 1:1000