

# CYP27A1 antibody



Catalog Number: 109689

## Product name

CYP27A1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200


## Immunogen

[Human CYP27A1 Recombinant protein \(His tag\)](#)

## Validations

## Specificity

Human; other species not tested.




## Antibody description

CYP27A1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HepG2 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 60 kDa

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:109689(CYP27A1 antibody) at dilution of 1:1000

## Preparation

This antibody was obtained by immunization of CYP27A1 recombinant protein (Accession Number: NM\_000784). Purification method: Antigen affinity purified.



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded human kidney using Catalog No:109689(CYP27A1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109689(CYP27A1

# CYP27A1 antibody



Catalog Number: 109689

---

antibody) at dilution of 1:100 (under 40x lens)