

# BCL11A antibody

Catalog Number: 108389

---

## Product name

BCL11A antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

BCL11A Rabbit Polyclonal antibody. Positive IHC detected in human gliomas.

## Preparation

This antibody was obtained by immunization of BCL11A recombinant protein (Accession Number: NM\_018014). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108389(BCL11A antibody) at dilution of 1:50 (under 10x lens)