

HSD17B13 Polyclonal Antibody



Catalog Number: 164789

Product name

HSD17B13 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to HSD17B13

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human HSD17B13 (NP_835236.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - HSD17B13 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HSD17B13 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure




Immunohistochemistry - HSD17B13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using HSD17B13 Antibody at dilution of 1:100 (40x lens).

HSD17B13 Polyclonal Antibody




Catalog Number: 164789



Immunohistochemistry - HSD17B13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using HSD17B13 Antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - HSD17B13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using HSD17B13 Antibody at dilution of 1:100 (40x lens).