

Anti-Cytochrome C Oxidase subunit Vic Rabbit antibody



Catalog Number: 166789

Product name

Anti-Cytochrome C Oxidase subunit Vic Rabbit antibody

Dilutions

WB: 1/1000

IHC: 1/20

Specificity

Human, Mouse, Rat

Antibody description

Cytochrome C Oxidase subunit Vic Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Cytochrome C Oxidase subunit Vic

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P

Validations



Immunohistochemistry of Cytochrome C Oxidase subunit Vic in paraffin-embedded Human colon cancer tissue using Cytochrome C Oxidase subunit Vic Rabbit pAb at dilution 1/20



Western blot detection of Cytochrome C Oxidase subunit Vic in 3T3, HeLa cell lysates using Cytochrome C Oxidase subunit Vic Rabbit pAb(1:1000 diluted). Predicted band size: 9kDa. Observed band size: 12kDa.