

VPS45 antibody



Catalog Number: 116789

Product name

VPS45 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS45 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue, human heart tissue, human testis tissue. Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of VPS45 recombinant protein (Accession Number: NM_007259). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116789(VPS45 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of HepG2 cells, using VPS45 antibody Catalog No:116789 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).