

VNN1 Polyclonal Antibody



Catalog Number: 164099

Product name

Rabbit IgG

VNN1 Polyclonal Antibody

Applications

Specificity

WB, IF

Human, Mouse

Dilutions

Antibody description

WB 1:500 - 1:2000

IF 1:50 - 1:100

Polyclonal antibody to VNN1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 300-415 of human VNN1 (NP_004657.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations



Western blot - VNN1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using VNN1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.