

Seladin 1 antibody



Catalog Number: 115099

Product name

Seladin 1 antibody

Immunogen

[Human DHCR24 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Seladin 1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HepG2 cells, MCF7 cells, mouse lung tissue. Positive IHC detected in human breast cancer tissue, human liver tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of Seladin 1 recombinant protein (Accession Number: NM_014762). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115099(DHCR24 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115099(DHCR24 antibody) at dilution of 1:100 (under 25x lens)