

Product name

CSN3 antibody

Specificity

Human; other species not tested.

Antibody description

CSN3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue.

Preparation

This antibody was obtained by immunization of CSN3 recombinant protein (Accession Number: XM_011531620). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:10-1:100

Validations

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109599(CSN3 Antibody) at dilution of 1:25 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109599(CSN3 Antibody) at dilution of 1:25 (under 40x lens)