

PHACTR4 antibody



Catalog Number: 113799

Product name

PHACTR4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PHACTR4 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue, A431 cells, human liver tissue. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 71-85 kDa

Preparation

This antibody was obtained by immunization of PHACTR4 recombinant protein (Accession Number: XM_017002168). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:113799(PHACTR4 antibody) at dilution of 1:400



Immunofluorescent analysis of A431 cells using Catalog No:113799(PHACTR4 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)