

Anti-CDK2 Rabbit antibody



Catalog Number: 169002

Product name

Anti-CDK2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CDK2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Cdk2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

1:2000


1:20

1:20


1:20

1:20

Validations



Immunocytochemistry of CDK2(green) in HeLa cells using CDK2 Rabbit mAb at dilution 1/200, and DAPI(blue)



and DAPI(blue)

Western blot detection of CDK2 in Jurkat, C6, CHO-K1, HeLa cell lysates using CDK2 Rabbit mAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.