

INMT antibody



Catalog Number: 111802

Product name

INMT antibody

Immunogen

[Human INMT Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

INMT Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue. Observed molecular weight by Western-blot: 29kd

Preparation

This antibody was obtained by immunization of INMT recombinant protein (Accession Number: NM_006774). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:111802(INMT antibody) at dilution of 1:300