

Syntenin 2 antibody



Catalog Number: 115802

Product name

Syntenin 2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Syntenin 2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human gliomas tissue. Positive WB detected in mouse brain tissue, human liver tissue, mouse heart tissue. Observed molecular weight by Western-blot: 35-37kd

Preparation

This antibody was obtained by immunization of Syntenin 2 recombinant protein (Accession Number: NM_015685). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115802(SDCBP2 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:115802(SDCBP2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog

Syntenin 2 antibody



Catalog Number: 115802

No:115802(SDCBP2 Antibody) at dilution of 1:50 (under 40x lens)