

SCO1 antibody



Catalog Number: 115012

Product name

SCO1 antibody

Specificity

Human; other species not tested.

Antibody description

SCO1 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of SCO1 recombinant protein (Accession Number: NM_004589). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115012(SCO1 antibody) at dilution of 1:400