

MALT1 antibody



Catalog Number: 107412

Product name

MALT1 antibody

Immunogen

[Human MALT1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

MALT1 Mouse Monoclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of MALT1 recombinant protein (Accession Number: NM_173844). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with (MALT1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:107412(MALT1 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)