

Anti-JNK1 Rabbit antibody



Catalog Number: 168812

Product name

Anti-JNK1 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

JNK1 Rabbit monoclonal antibody

Dilutions

1:1000

Preparation

Antigen: Recombinant protein of human JNK1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of JNK1 in Jurkat,Rat Brain,C6 cell lysates using JNK1 Rabbit mAb(1:1000 diluted).Predicted band size:48kDa.Observed band size:46kDa.