

C2orf3 antibody



Catalog Number: 108712

Product name

C2orf3 antibody

Immunogen

[Human C2orf3 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C2orf3 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 89 kDa

Preparation

This antibody was obtained by immunization of C2orf3 recombinant protein (Accession Number: NM_003203). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108712(C2orf3 Antibody) at dilution of 1:600