

CDC25A Polyclonal Antibody



Catalog Number: 160822

Product name

CDC25A Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CDC25A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 235-484 of human CDC25A (NP_963861.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - CDC25A Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat liver using CDC25A antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CDC25A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using CDC25A antibody at dilution of 1:100



(40x lens).


Immunohistochemistry - CDC25A Polyclonal Antibody

CDC25A Polyclonal Antibody



Catalog Number: 160822

Immunohistochemistry of paraffin-embedded rat heart using CDC25A antibody at dilution of 1:100




(40x lens).

Immunohistochemistry - CDC25A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using CDC25A antibody at

dilution of 1:100 (40x lens).



Immunohistochemistry - CDC25A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using CDC25A antibody at dilution of 1:100 (40x lens).