

GPR19 antibody



Catalog Number: 111122

Product name

GPR19 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GPR19 Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue.

Preparation

This antibody was obtained by immunization of GPR19 recombinant protein (Accession Number: XM_047428742). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111122(GPR19 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:111122(GPR19 Antibody) at dilution of 1:100 (under 40x lens).