

OVOL1 antibody



Catalog Number: 113522

Product name

OVOL1 antibody

Immunogen

[Human OVOL1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

OVOL1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse pancreas tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of OVOL1 recombinant protein (Accession Number: NM_004561). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113522(OVOL1 antibody) at dilution of 1:400



IP Result of anti-OVOL1 (IP:Catalog No:113522, 4ug; Detection:Catalog No:113522 1:500) with mouse brain tissue lysate 4000ug.