

DHRS11 antibody

Catalog Number: 109922

Product name

DHRS11 antibody

Immunogen

[Human DHRS11 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DHRS11 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, COLO 320 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 25-28kd

Preparation

This antibody was obtained by immunization of DHRS11 recombinant protein (Accession Number: NM_024308). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109922(DHRS11 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:109922(DHRS11 Antibody) at dilution of 1:50(under 40x lens)