

Follistatin antibody



Catalog Number: 107222

Product name

Follistatin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Follistatin Mouse Monoclonal antibody. Positive IF detected in Hela cells. Positive WB detected in Hepg2 cells, HeLa cells, Raji cells. Observed molecular weight by Western-blot: 32 / 35/ 40 kDa

Preparation

This antibody was obtained by immunization of Follistatin recombinant protein (Accession Number: NM_013409). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107222(FST antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using FST antibody Catalog No: at 1:25 dilution and Rhodamine-labeled goat anti-mouse IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).