

SAP18 Polyclonal Antibody



Catalog Number: 164322

Product name

SAP18 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SAP18

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 21-172 of human SAP18 (NP_005861.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:200 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - SAP18 Polyclonal Antibody


Western blot analysis of extracts of SW480 cells, using SAP18 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - SAP18 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using SAP18 antibody at dilution of 1:100



(40x lens).

SAP18 Polyclonal Antibody



Catalog Number: 164322

Immunohistochemistry - SAP18 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using SAP18 antibody at dilution of 1:100 (40x lens).