

ACADSB antibody



Catalog Number: 107722

Product name

ACADSB antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ACADSB Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human breast cancer. Observed molecular weight by Western-blot: 45-47kd

Preparation

This antibody was obtained by immunization of ACADSB recombinant protein (Accession Number: NM_001609). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107722(ACADSB antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107722(ACADSB antibody) at dilution of 1:100 (under 10x lens)