

ANGPTL3 Polyclonal Antibody



Catalog Number: 160132

Product name

ANGPTL3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ANGPTL3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 231-460 of human ANGPTL3 (NP_055310.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:100 - 1:200

Validations



Western blot - ANGPTL3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ANGPTL3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - ANGPTL3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using ANGPTL3 Antibody at dilution of 1:100 (40x lens).