

Anti-PKM2 Rabbit antibody



Catalog Number: 167132

Product name

Anti-PKM2 Rabbit antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

PKM2 Rabbit polyclonal antibody

Preparation

Antigen: CVTEVENGGSLGSKKGVNLP

Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of PKM2 in COS7,C6,3T3 and HeLa cell lysates using PKM2 Rabbit pAb (1:1000 diluted).Predicted band size:60kDa.Observed band size:60kDa.