

RALYL antibody



Catalog Number: 114532

Product name

RALYL antibody

Immunogen

[Human RALYL Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

RALYL Rabbit Polyclonal antibody. Positive IHC detected in mouse brain tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of RALYL recombinant protein (Accession Number: XM_017013065). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114532(RALYL antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114532(RALYL Antibody) at dilution of 1:50 (under 40x lens)